

## Superior Closure Prevents Fluid Exposure

- ▶ Provides Leak Tight Seal
- ▶ Offers One-Handed Manipulation
- ▶ Accommodates Endourological  
0-5F Instruments

# Sureseal<sup>®</sup> II

Endoscopic Valve



# Sureseal® II Endoscopic Valve



## Superior Closure

*Keeps urologist dry. Works with a variety of scopes and the Applied Forté ureteral access sheaths. Accommodates instruments from 0-5F.*

## Enhanced Handling

*Allows passage of instruments into working port without additional manipulation of valve. Enables one-handed manipulation of scope.*

## Latex Free

*Protects latex sensitive patients.*

## The Seal That Keeps You Dry

The Sureseal II endoscopic valve combines maximum irrigation capability with a leak-tight seal around 0-5F stone baskets, guidewires, laser fibers and other endourological instruments.

REORDER NO.	VALVE OPENING	PRESSURE RATING	IRRIGATION PORT	SCOPE PORT	QUANTITY
Seal-005	Passive, continuously reliable bidirectional from 0-5F (.065")	20psi or 1034mm Hg	Luer-lock type fitting	Standard luer-lock with ball adapter	10/BOX (includes valve and ball adapter)



© 2013 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a \* are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01077G

<http://contact.appliedmedical.com/>

For more information [www.appliedmedical.com/accessories](http://www.appliedmedical.com/accessories)

