



# Inzii<sup>®</sup>

## Retrieval Systems

A large teal graphic element consisting of a diagonal line that starts from the top left and extends towards the bottom right, creating a triangular shape on the left side of the page.

# COMPLETE SPECIMEN RETRIEVAL BAG OFFERING

Exceptional durability ensures specimen security and minimizes risk of contamination

Innovative design allows introduction through smaller ports, advancing minimally invasive procedures



### **Inzii 12/15mm**

15mm bag in a 12mm\* device



### **Inzii 10mm**

Standard specimen retrieval bag

## Procedural Applications


Facilitates a wide range of laparoscopic procedures, including:



### LAPAROSCOPIC

Cholecystectomy  
Oophorectomy  
Appendectomy

Nephrectomy  
Splenectomy  
Sleeve Gastrectomy



## Guide Bead

Allows reopening and closure of bag, facilitating multiple specimen retrieval



## Inzii Universal 10mm bag in a 5mm\* device

Reinforced tip design ensures specimen security and minimizes risk of contamination

Bag remains open and intact throughout specimen capture





## Inzii Universal Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD003	5mm*	180ml	10/Box



## Inzii 10mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD001	10mm	225ml	10/Box



## Inzii 12/15mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD004	12mm*	1600ml	5/Box

\* Applied Medical Kii trocar

Visit [www.appliedmedical.com/inzii](http://www.appliedmedical.com/inzii)  
for more information

© 2017 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design and marks designated with a  are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States and/or the European Union. SC01638G

