



Applied Scope Warmer

A large teal graphic element consisting of a diagonal line that divides the page into two sections. The upper section is white and contains the main title and two lines of text. The lower section is a solid teal color.

A SIMPLE SOLUTION FOR BETTER VISIBILITY

Simple to use
Minimizes lens fogging

The Answer to Laparoscopic Fogging

Eliminates the need for chemical gels or wipes, retains heat throughout lengthy procedures and minimizes the need for repeated scope cleansing and defogging.

Scope Warmer Seal

Retains heat throughout lengthy procedures

Scope Warmer Bottle

Minimizes fogging by pre-warming the laparoscope to body temperature

Scope Warmer Base

Stabilizes the system when a scope is inserted





Applied Scope Warmer

Reorder No.	Description	Qty/Box
C3001	Scope Warmer Bottle* Reusable Volume Capacity = 500mL	1 EACH
C3002	Scope Warmer Base* Reusable	1 EACH
C3101	Scope Warmer Seals Disposable, Single Use Individually packaged in sterile pouches	20/BOX

Recommended Sterilization Method: Autoclave or STERRAD **100S.
** STERRAD is a registered trademark of Advanced Sterilization Products.
Please see Instructions for Use (IFU) for complete information.
Scope is not included.

Visit www.appliedmedical.com/laparoscopic
for more information.

© 2017 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks
of Applied Medical Resources Corporation, registered in one or more of the following countries:
Australia, Canada, Japan, South Korea, the United States, and/or the European Union.
SC01094J

