

A Simple Solution for Better Visibility

- ▶ Simple to use
- ▶ Minimizes lens fogging

Applied Scope Warmer



Applied
Medical
SCOPE WARMER
REF C3001

CE

Applied
Medical 
A New Generation Medical Device Company

Applied Scope Warmer



Scope Warmer Seal

Retains heat throughout lengthy procedures

Scope Warmer Bottle

Minimizes fogging by pre-warming the laparoscope to body temperature

Scope Warmer Base

Stabilizes the system when a scope is inserted

The Solution to Laparoscopic Fogging

Eliminates the need for chemical gels or wipes, retains heat throughout lengthy procedures, and minimizes the need for repeated scope cleansing and defogging.

| REORDER NO. | DESCRIPTION | QUANTITY |
|-------------|--|----------|
| C3001 | Scope Warmer Bottle* Reusable Volume Capacity = 500mL | 1 EACH |
| C3002 | Scope Warmer Base* Reusable | 1 EACH |
| C3101 | Scope Warmer Seals Disposable, Single Use Individually packaged in sterile pouches | 20/BOX |

* Recommended Sterilization Method: Autoclave or STERRAD**1100S.

** STERRAD is a registered trademark of Advanced Sterilization Products.
Please see Instructions for Use (IFU) for complete information.
Scope is not included.



© 2013 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a * are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01094H

<http://contact.appliedmedical.com/>

For more information www.appliedmedical.com/accessories

