




AeroStat[®]

Cholangiography Catheter

ONE HAND ONE STEP

Easy duct access
Self-anchoring
One-handed activation



Easy Duct Access

Retention cone eliminates problems traditionally associated with laparoscopic cholangiography, such as clipping, balloon purging and potential breakage.

Self-Anchoring

Soft 6F-20F expandable cone and 3F tapered rigid tip atraumatically engages and occludes the cystic duct, eliminating the need for clips or graspers.

One-Handed Activation

Ratcheting handle, with one-handed activation, makes placement fast and easy.





AeroStat Cholangiography Catheter

Reorder No.	Description	Shaft Diameter	Length
C1002	Rigid Tip with Percutaneous Introducer	6F	33cm

Visit www.appliedmedical.com/laparoscopic
for more information.