

**One Hand. One Step.**

- ▶ Easy duct access
- ▶ Self-anchoring
- ▶ One-handed activation

# AeroStat<sup>®</sup>

Cholangiography Catheter



**Applied**  
Medical




*A New Generation Medical Device Company*

# AeroStat® Cholangiography Catheter

## Easy Duct Access

*Retention cone eliminates problems traditionally associated with laparoscopic cholangiography, such as clipping, balloon purging and potential breakage.*




## Self-anchoring

*Soft 6F-20F expandable cone and 3F tapered rigid tip atraumatically engages and occludes the cystic duct, eliminating the need for clips or graspers.*

## One-handed Activation

*Ratcheting handle, with one-handed activation, makes placement fast and easy.*



REORDER NO.	DESCRIPTION	SHAFT DIAMETER	LENGTH
C1002	Rigid Tip with Percutaneous Introducer	6F	33cm



© 2015 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a \* are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01096F

<http://contact.appliedmedical.com/>

For more information [www.appliedmedical.com/accessories](http://www.appliedmedical.com/accessories)

**Applied**  
Medical

PRO TERRA