

# STERILE SPECIMEN COLLECTION SWAB

6-Inch Lattice-Tip Nasopharyngeal Swab  
with 70mm Breakpoint



As we navigate through these unprecedented times, Applied Medical is committed to rapidly responding to the needs of the healthcare community. We are all in this together, fighting against COVID-19.

# 6-INCH LATTICE-TIP NASOPHARYNGEAL SWAB WITH 70mm BREAKPOINT

## PRODUCT DESCRIPTION

Applied Medical's lattice-tip sterile specimen collection swab is indicated for the collection and transportation of nasopharyngeal clinical samples to support COVID-19 testing. The swab consists of a handle and a lattice-tip cover that is securely bonded to a medical-grade polystyrene microbrush-shaped tip, to support enhanced absorption.

## FEATURES

- Sterile, individually packaged swabs
- Made in the USA with materials sourced domestically

## DIMENSIONS

Overall Length: 5.9 inches

### Tip

- Design: Microbrush with lattice cover
- Material: Polystyrene and polyester with adhesive
- Length: 0.757 inches
- Diameter: 0.014 inches

### Handle

- Material: Polystyrene
- Length: 3.096 inches
- Diameter: 0.098 inches

Breakpoint: 70mm

## ORDER INFORMATION

Model #	Description	Qty/Box
NPS20	Sterile specimen collection swab	200*

\*Controlled by weight

Applied Medical is pleased to offer nasopharyngeal swabs direct to U.S. healthcare facilities and government agencies. For more information or to inquire about placing an order, please contact [COVIDResponse@appliedmedical.com](mailto:COVIDResponse@appliedmedical.com).

