

FACE SHIELD

FOR MEDICAL USE



As we navigate through these unprecedented times, Applied Medical is committed to rapidly responding to the needs of the healthcare community. We are all in this together to fight against COVID-19.


FACE SHIELD

PRODUCT DESCRIPTION

The Applied Medical face shield was designed and developed to meet the needs of our healthcare partners on the front lines of this crisis. It is carefully crafted with materials on hand to help protect those of you serving others. The face shields may vary in colors or materials as we focus on the functionality, comfort and quality of the product.

FEATURES

- Foam headband in an assortment of colors (white, yellow and charcoal)
- Elastic strap to accommodate a variety of head sizes
- Clear shield for optimal view
- Industrial staples secure the strap to the shield
- Fits over goggles or glasses
- Not made with natural rubber latex
- Made in the USA with materials sourced domestically
- Not recommended for use in the presence of high intensity heat sources or flammable gases
- Disposable
- Non-sterile
- Single use



BODY CONTACTING MATERIALS

- Foam
- Polycarbonate or polyethylene (PETG)
- Dynaflex
- Steel

ORDER INFORMATION

Model #	Description	Face Shield Length	Qty
FS001	Face shield for medical use, one size fits most	9"	100/Box (packaged 10/pack)

Applied Medical is pleased to offer face shields direct to US healthcare facilities and government agencies. For more information or to inquire about placing an order, please contact COVIDResponse@appliedmedical.com.



This face shield is not intended to prevent specific diseases or infections. This product has not been FDA cleared or approved. It has been authorized by FDA under an EUA for use by healthcare providers as personal protective equipment. This face shield is only authorized for the duration of the declaration that circumstances justifying the authorization of emergency use under Section 546(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.

© 2020 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union.
297966-EN-USA-B

