



# Epix<sup>®</sup> Universal

Clip Applier

# ONE SOLUTION FOR LIGATION



## Low Profile Jaw

Allows for maximum visualization and clearance around the tissue to be ligated

Allows for greater clip retention through the application of a true M/L clip



## 34cm Shaft

Maximizes instrument reach

## Advanced Clip Delivery System

Offers reliable, smooth device actuation

## Ergonomic Handle Design


Provides optimal control during clip loading and closure



True M/L clip

## Superior Clip Retention

Features TiGold® clip technology, ensuring reliable closure and secure retention





## Epix Universal Clip Applier

Reorder No.

Qty/Box

*Epix Universal (contains 20 M/L TiGold clips)*

CA500\*

3/Box

\* Epix Universal clip applier is compatible with Applied Medical Kii® 5mm and larger trocars.

Visit [www.appliedmedical.com/epix](http://www.appliedmedical.com/epix)  
for more information

Devices listed may not be approved in all markets.  
Please contact your field implementation team member for more information on availability.

© 2019 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. All of the other trademarks used herein are the property of their respective owners. 194619-EN-USA-A

