



Forté[®]

Ureteral Access Sheaths

Forté

Ureteral Access Sheaths

Advanced Technology, Clinically Proven

Forté ureteral access sheaths provide a continuous working channel, simplifying ureteroscopic procedures.

Forté sheaths have been clinically proven to protect the urinary system, reduce operating room time, eliminate the need for balloon dilators, and aid in washing out stone fragments.^{1, 2, 3}



Seamless Transition

The sheath-and-dilator combination facilitates placement.

Continual Access

Providing continual access through a sheath reduces procedural time.
Using a sheath protects the scope and instruments, reducing the need for costly repairs.

Less Trauma

The risk of trauma decreases when the ureter is protected by a sheath during ureteroscopy.
The hydrophilic coating facilitates advancement.

Gentle Dilation

The tapered dilator gently traverses the ureteral orifice and replaces expensive balloon dilators.
The dilator ensures easy insertion and minimizes trauma.



Forté Ureteral Access Sheaths

Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7266	10F	12F–16F	20cm	6
B7276	10F	12F–16F	28cm	6
B7216	10F	12F–16F	35cm	6
B7296	10F	12F–16F	45cm	6
B7286	10F	12F–16F	55cm	6



Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7066	12F	14F–18F	20cm	6
B7076	12F	14F–18F	28cm	6
B7016	12F	14F–18F	35cm	6
B7096	12F	14F–18F	45cm	6
B7086	12F	14F–18F	55cm	6



Reorder No.	Dilator	Sheath	Length	Qty/Box
<i>AxP</i>				
B7166	14F	16F–18F	20cm	6
B7176	14F	16F–18F	28cm	6
B7116	14F	16F–18F	35cm	6

1. Auge BK, Pietrow PK, Lallas CD, Raj GV, Santa-Cruz RW, Preminger GM. Ureteral access sheath provides protection against elevated renal pressures during routine flexible ureteroscopic stone manipulation. *J Endourol.* 2004;18(1):33-36. doi:10.1089/089277904322836631.

2. Kourambas J, Byrne RR, Preminger GM. Does a ureteral access sheath facilitate ureteroscopy?. *J Urol.* 2001;165(3):789-793. doi:10.1016/S0022-5347(05)66527-5.

3. Ajayi L, Lillas L, Montgomery B, Naerger H, Palfrey E. The use of ureteral access sheaths and flexible ureteroscopy; a marriage in heaven? *J Endourol.* 2001;15(suppl 1).

Visit [appliedmedical.com/forte](https://www.appliedmedical.com/forte) for more information.

© 2022 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.