



Inzii[®]

Retrieval Systems

A large teal graphic element consisting of a diagonal line that starts from the top left and extends towards the bottom right, creating a triangular shape on the left side of the page.

COMPLETE SPECIMEN RETRIEVAL BAG OFFERING

Reinforced tip design ensures specimen security and
minimizes risk of contamination

Innovative design allows introduction through smaller
ports, advancing minimally invasive procedures

Bag remains open and intact throughout
specimen capture

Inzii Universal 10mm bag in a 5mm* device

Minimizes the need for large ports and enhances versatility

Extra-large angled design eases specimen delivery into bag





Inzii 12/15mm

15mm bag in a 12mm* device



Inzii 10mm

Standard specimen retrieval bag

Procedural Applications


Facilitates a wide range of laparoscopic procedures, including:



LAPAROSCOPIC

Cholecystectomy
Oophorectomy
Appendectomy

Nephrectomy
Splenectomy
Sleeve Gastrectomy



Guide Bead

Allows reopening and closure of bag, facilitating multiple specimen retrieval





Inzii Universal Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD003	5mm*	180ml	10/Box



Inzii 10mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD001	10mm	225ml	10/Box



Inzii 12/15mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD004	12mm*	1600ml	5/Box

* Applied Medical Kii trocar

Visit www.appliedmedical.com/inzii
for more information

© 2017 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States and/or the European Union. SC01638H

